



CO (ppm)

%COHB

20 & Above	Heavy Smoker	Heavy Smoker	3.20+
20			3.20
19			3.04
18			2.88
17			2.72
16	Smoker	Smoker	2.56
15			2.40
14			2.24
13			2.08
12			1.92
11			1.76
10			1.60
9	Light Smoker	Light Smoker	1.44
8			1.28
7			1.12
6			0.96
5			0.80
4			0.64
3	Non-Smoker	Non-Smoker	0.48
2			0.32
1			0.16